

DIGITAL PERSONAL SCALE

7810



USER MANUAL

Read the user manual carefully before use.

1. FEATURES:

Maximum Capacity:	150 kg / 330 lb
Divison:	0.1 kg / 0.2 lb
Display:	LCD-display, 4 digits
Battery:	1x CR2032
Working temperature:	0° - 35° C
Units:	kg / lb / st
Auto-Off-Function:	after 8 seconds of non-use

2. OPERATION INSTRUCTIONS

Please remove the plastic foil from the battery before use.

Place the scale on an even floor. An uneven floor will affect the accuracy of measurement.
Stand on the scale in balance and start weighing. The measurement will appear for 8 seconds.
Step off the scale. The scale will turn off automatically.
Replace the battery if the display shows "LO" or "- - -".
Change units by switching the button on the bottom of the scale.
Keep the scale dry.

3. ATTENTION

Please do not load the scale after use in order to avoid unwanted unloading of battery.
Try to stand in the middle of the platform to avoid slipping off.
Stand still on the scale to avoid inaccurate reading.
Keep the scale dry. Do not clean the scale under running water.
Remove battery when unused for a longer period.

4. ERRORS & REMEDIES

Error code	Reason	Remedies
"LO" or "- - -"	battery is empty	change battery
"ERR"	scale is unstable	place the scale on an even floor
"_ _ _"	if appears during weighing, max. capacity has been exceeded	step off the scale

Soehnle Industrial Solutions GmbH
Gaildorfer Straße 6
DE-71522 Backnang
GERMANY

Telefon: +49 (0) 71 91 / 34 53 - 220
Telefax: +49 (0) 71 91 / 34 53 - 211
info@soehnle-professional.com
www.soehnle-professional.com