

# **RD-953** innerScan.

### **BLUETOOTH BODY COMPOSITION MONITOR WITH EXCLUSIVE MULTI FREQUENCY BIA TECHNOLOGY**



FDA Cleared "Multi Frequency" BIA Technology

Integrated Bluetooth 4.0 long battery life allowing instant data transmission to your iPhone 4S and above or Android 4.4 and above

Easy-to-read measurements with instant color-coded assessments

User recognition - just step on and the scale will start taking a measurement

11 Full-Body Measurements, including muscle quality score

Graph



### Free Tanita Healthy Edge® Mobile App

Automatically upload your measurements to Tanita Healthy Edge® Mobile App.


 $\checkmark$ iPhone 4S or above  $\checkmark$ Android 4.4 or above

Provides a powerful tool to track your progress with easy-to-read graphs and helpful assessments. Set your own targets and achieve your personal fitness and wellness goals.

## **RD-953** innerScand BLUETOOTH BODY COMPOSITION MONITOR

Tanita, the designer, manufacturer & distributor of the world's most accurate, radio wireless body composition scales. Since creating the world's first stand-on body composition monitor in 1992, Tanita Bioelectrical Impedance Analysis (BIA) technology has been used by health and well-being experts around the world.

Tanita has produced the highest standard in fitness and wellness assessment. The RD-953 incorporates the very best *FDA Cleared*, medical grade "*Multi Frequency*", BIA technology, allowing you to have the most accurate body composition accuracy in the comfort of your home.




#### What is "Multi Frequency"?

Research has shown that using multiple bio-electrical impedance (BIA) frequencies provides essential data of a person's *intracellular* and *extracellular status*. This advanced technology allows greater accuracy when calculating body composition measurements, including muscle quality score.

#### YOUR MEASUREMENTS ARE ASSESSED AND SHOWN IN 3 SIMPLE LEVELS





### Muscle Quality Score




Higher density of muscle fiber with little fat, water and connective tissue, etc.

⇒ Maintain regular exercise to ensure your Muscle Quality stays within the healthy range.



#### **Specifications:**

- ✓ Bluetooth 4.0 Low Energy
- ✓ Capacity: 440lb (200 kg)
- ✓ Weight increments: 0-200lb / 0.1lb, 200-440lb / 0.2lb 0-100kg / 0.05kg, 100-200kg / 0.1kg
- ✓ Body Fat Increments: 0.1%
- ✓ Tempered glass platform with transparent electrodes
- ✓ 4 user memory each with it's own color display
- ✓ Color Coded Indicators
- ✓ Weight only function
- ✓ Guest Mode & Athlete Mode
- ✓ Power Supply: 4 x AA (included)
- ✓ Dimensions: 11.7" x 12.9" x 1.3"
- ✓ Weight: 4.6lb





Muscle fiber is thinner and fewer as compared with other substances

⇒ Increase exercise on a regular basis to improve your Muscle Quality Score



2625 S. Clearbrook Drive / Arlington Heights, IL 60005 / Phone: 847-640-9241 / Fax: 847-640-7978 / E-mail: 4health@tanita.com / Web: www.Tanita.com